

## pH Tester with Replaceable Electrode



At the cost of just a few rolls of litmus paper, you can have the industry's best pocket-sized pH tester.

The *Checker*<sup>®</sup> supplies you with fast and accurate readings from 0 to 14 pH with a resolution of 0.01 pH.

This compact tester features a large easy-to-read LCD and simple to perform 2-point calibration.

The pH electrode of the *Checker*<sup>®</sup> can be replaced. You only need to unscrew the electrode from the meter body and screw on a new one.

The *Checker*<sup>®</sup> is fast, accurate, light-weight and with 3000 hours of battery life, you don't have to worry about frequent battery replacement.

Testing pH has never been easier!

### Specifications

HI 98103 <i>Checker</i> <sup>®</sup>	
Range	0.00 to 14.00 pH
Resolution	0.01 pH
Accuracy (@20°C/68°F)	±0.2 pH
Calibration	manual, 2 points
Electrode	HI 1270 (included)
Battery Type / Life	2 x 1.5V / approx. 3000 hours of continuous use
Environment	0 to 50°C (32 to 122°F); RH max 95%
Dimensions	66 x 50 x 25 mm (2.6 x 2.0 x 1.0") - without probe
Weight	50 g (1.8 oz.)

### Ordering Information

HI 98103 (*Checker*<sup>®</sup>) is supplied complete with HI 1270 pH electrode, batteries and instructions.

### Accessories

HI 1270	pH electrode with screw-type connector	HI 7004M	pH 4.01 buffer solution, 230 mL bottle
HI 77400P	pH 4 and pH 7 buffer sachets (20 mL, 5 pcs./each)	HI 7007M	pH 7.01 buffer solution, 230 mL bottle
HI 70004P	pH 4.01 buffer solution, 20 mL sachet, 25 pcs.	HI 7010M	pH 10.01 buffer solution, 230 mL bottle
HI 70007P	pH 7.01 buffer solution, 20 mL sachet, 25 pcs.	HI 70300M	Electrode storage solution, 230 mL bottle
HI 70010P	pH 10.01 buffer solution, 20 mL sachet, 25 pcs.	HI 7061M	Cleaning solution, 230 mL bottle

For a complete range of calibration, cleaning and maintenance solutions, see section F.