

1413 $\mu\text{S}/\text{cm}$ Calibration Solution

The 1413 $\mu\text{S}/\text{cm}$ solution is the one best suited for general use.

This solution is also available in combined kits with the pH 7 buffer for the easy calibration of multiparameter instruments.

It is also available in different sized bottles and in single-dose, ready-to-use sachets.

The HI 8031L solution is provided in an opaque bottle according to FDA (Food & Drug Administration) regulations, which prevents the reagent from damage due to extended exposure to light.

Our wide range of calibration solutions also includes solutions provided with a certificate of analysis, to satisfy the requirements of any application from the farm to the factory.



HI 6031: High Precision Standard Solution

- **Standard Solution at 1413 $\mu\text{S}/\text{cm}$ (@25°C)**



Solution Bottle

Code	EC Value @25°C	Size	Package	FDA Bottle	Certificate of Analysis
HI 6031	1413 $\mu\text{S}/\text{cm}$ ($\mu\text{mho}/\text{cm}$)	500 mL	1 bottle		•
HI 7031/1G	1413 $\mu\text{S}/\text{cm}$ ($\mu\text{mho}/\text{cm}$)	1 Gallon (3.78 L)	1 bottle		
HI 7031L	1413 $\mu\text{S}/\text{cm}$ ($\mu\text{mho}/\text{cm}$)	500 mL	1 bottle		
HI 7031L/C	1413 $\mu\text{S}/\text{cm}$ ($\mu\text{mho}/\text{cm}$)	500 mL	1 bottle		•
HI 7031M	1413 $\mu\text{S}/\text{cm}$ ($\mu\text{mho}/\text{cm}$)	230 mL	1 bottle		
HI 8031L	1413 $\mu\text{S}/\text{cm}$ ($\mu\text{mho}/\text{cm}$)	500 mL	1 bottle	•	

Solution Sachet

Code	EC Value @25°C	Size	Package	Certificate of Analysis
HI 70031C	1413 $\mu\text{S}/\text{cm}$ ($\mu\text{mho}/\text{cm}$)	20 mL	25 sachets	•
HI 70031P	1413 $\mu\text{S}/\text{cm}$ ($\mu\text{mho}/\text{cm}$)	20 mL	25 sachets	
HI 77100C	1413 $\mu\text{S}/\text{cm}$ e pH 7.01	20 mL	20 sachets (10 each)	•
HI 77100P	1413 $\mu\text{S}/\text{cm}$ e pH 7.01	20 mL	20 sachets (10 each)	