

## Accessories

### Beaker

Rugged plastic beakers available in 2 different sizes. For use during testing with HANNA instruments® chemical test kits and photometers. Anti-spill ported caps are available separately.

- HI 740036P 50 mL plastic beaker (10 pcs)
- HI 740037P 20 mL plastic beaker (10 pcs)
- HI 740032P Cap for 20 mL plastic beaker (10 pcs)
- HI 740034P Cap for 50 mL plastic beaker (10 pcs)

### Blades for Meat pH Electrodes

These stainless steel blades, in conjunction with FC 230, FC 231 and FC 232 pH electrodes, allow easy meat penetration. Two models are available for different depths.

- FC 098 Stainless steel blade for meat pH electrode, 20 mm (0.8")
- FC 099 Stainless steel blade for meat pH electrode, 35 mm (1.4")

### Boot, Shockproof

These special boots have been designed to protect your instruments and prevent damage due to accidental impact, particularly in field applications. They can withstand high temperature (up to 200°C / 392°F), ozone and environmental phenomena, and are available in 2 colors.

Shockproof rubber boot for portable meters with HI 991001 casing.

- HI 710007 Shockproof rubber boot, blue
- HI 710008 Shockproof rubber boot, orange

Shockproof rubber boot for portable meters with HI 93711 casing.

- HI 710009 Shockproof rubber boot, blue
- HI 710010 Shockproof rubber boot, orange

Shockproof silicone boot for portable meters with HI 8424NEW casing.

- HI 710015 Shockproof silicone boot, blue
- HI 710016 Shockproof silicone boot, orange

### Boron, Test Kit Reagents

see Section V

### Bottle for D.O. Analysis

Glass bottle, supplied with the HANNA instruments® chemical test kits and photometer for D.O. analysis. The bottle is provided with a glass stopper to ensure an air tight seal. Bottle and stopper can be also ordered separately.

- HI 740039 60 mL glass bottle for D.O. analysis
- HI 740041 60 mL glass bottle and stopper (6 pcs)
- HI 740042 60 mL glass bottle for D.O. analysis (6 pcs)
- HI 740043 Glass stopper for 60 mL glass bottle (6 pcs)

### Bottle, Dropper, Polyethylene

Provided with HANNA instruments® test kits where one-by-one drop dosing is required. Thanks to the removable plug, bottle refilling is very easy.

- HI 740060P 30 mL dropping bottle (10 pcs)

