

## Accessories



**HI 778P** Cable for electrodes with screw-cap connector, 100 m (333')

3 mm (0.12') dia. connection cable with a screw-cap connector at one end. The other end can be free or provided with a BNC, DIN or US standard connector.

- HI 7854/1 Connection cable with screw-cap connector, 1 m (3.3')
- HI 7854/3 Connection cable with screw-cap connector, 3 m (9.9')
- HI 7854/5 Connection cable with screw-cap connector, 5 m (16.5')
- HI 7854/10 Connection cable with screw-cap connector, 10 m (33')
- HI 7854/15 Connection cable with screw-cap connector, 15 m (49.5')
- HI 7855/1 Connection cable with screw & BNC connectors, 1 m (3.3')
- HI 7855/3 Connection cable with screw & BNC connectors, 3 m (9.9')
- HI 7855/5 Connection cable with screw & BNC connectors, 5 m (16.5')
- HI 7855/10 Connection cable with screw & BNC connectors, 10 m (33')
- HI 7855/15 Connection cable with screw & BNC connectors, 15 m (49.5')
- HI 7856/1 Connection cable with screw & US connectors, 1 m (3.3')
- HI 7856/3 Connection cable with screw & US connectors, 3 m (9.9')
- HI 7856/5 Connection cable with screw & US connectors, 5 m (16.5')
- HI 7856/10 Connection cable with screw & US connectors, 10 m (33')
- HI 7856/15 Connection cable with screw & US connectors, 15 m (49.5')
- HI 7857/1 Connection cable with screw & DIN connectors, 1 m (3.3')
- HI 7857/3 Connection cable with screw & DIN connectors, 3 m (9.9')
- HI 7857/5 Connection cable with screw & DIN connectors, 5 m (16.5')
- HI 7857/10 Connection cable with screw & DIN connectors, 10 m (33')
- HI 7857/15 Connection cable with screw & DIN connectors, 15 m (49.5')

3 mm (0.12') dia. connection cable with a BNC connector at one end. The other end can be provided with a BNC, DIN, US standard connector or spade lug.

- HI 7858/1 Connection cable with BNC & BNC connectors, 1 m (3.3')
- HI 7858/3 Connection cable with BNC & BNC connectors, 3 m (9.9')
- HI 7858/5 Connection cable with BNC & BNC connectors, 5 m (16.5')
- HI 7858/10 Connection cable with BNC & BNC connectors, 10 m (33')
- HI 7858/15 Connection cable with BNC & BNC connectors, 15 m (49.5')
- HI 7859/1 Connection cable with BNC & DIN connectors, 1 m (3.3')
- HI 7859/3 Connection cable with BNC & DIN connectors, 3 m (9.9')
- HI 7859/5 Connection cable with BNC & DIN connectors, 5 m (16.5')
- HI 7859/10 Connection cable with BNC & DIN connectors, 10 m (33')
- HI 7859/15 Connection cable with BNC & DIN connectors, 15 m (49.5')
- HI 7860/1 Connection cable with BNC & US connectors, 1 m (3.3')
- HI 7860/3 Connection cable with BNC & US connectors, 3 m (9.9')
- HI 7860/5 Connection cable with BNC & US connectors, 5 m (16.5')
- HI 7860/10 Connection cable with BNC & US connectors, 10 m (33')
- HI 7860/15 Connection cable with BNC & US connectors, 15 m (49.5')
- HI 7861/1 Connection cable with BNC conn. & spade lug, 1 m (3.3')
- HI 7861/3 Connection cable with BNC conn. & spade lug, 3 m (9.9')
- HI 7861/5 Connection cable with BNC conn. & spade lug, 5 m (16.5')
- HI 7861/10 Connection cable with BNC conn. & spade lug, 10 m (33')
- HI 7861/15 Connection cable with BNC conn. & spade lug, 15 m (49.5')



### **Cable, Extension, for K-Type Thermocouple Probes**

This coiled cable can be used to extend the cable of all K-type thermocouple probes. It comes with male and female connectors at the two ends, respectively.

**HI 766EX** Extension cable for K-type thermocouple probes, 1 m (3.3')

### **Cable for PC Connection**

RS232 serial cables to connect a PC to any HANNA instruments® meter provided with RS232 serial port.

PC Connection cable for bench meters.

HI 920010 Serial cable for PC connection (9 to 9-pin)

PC connection cable for portable meters.

HI 920011 Serial cable for PC connection (5 to 9-pin)

### **Calcium and Magnesium, Colorimeter Reagents**

see Section V

### **Calcium, Test Kit Reagents**

see Section V

### **Calibration Solution**

see Section F

### **Cap, Beaker**

see Beaker

### **Cap, Bottle**

see Bottle

### **Cap, Cuvet**

To be used with cuvetts for photometers and turbidity meters. They are notched to fit the groove on the meter, thus no light can enter and interfere with measurement.

HI 731325 Cap for cuvet for turbidity meters and HI 937xx series (4 pcs)

HI 731325N Cap for cuvet for HI 832xx series (4 pcs)

HI 731335 Cap for cuvet for HI 957xx series (4 pcs)

### **Cap, Protective**

Prevents tester and probe damage due to accidental impacts.

HI 740162P Protective cap for *Checktemp*® (5 pcs)

HI 740164 Protective cap for HI 1083

HI 740165 Protective cap for HI 1330

HI 740166 Protective cap for HI 1331

HI 740167 Protective cap for HI 2031

HI 740169 Protective cap for *Hydrocheck*®

HI 740170 Protective cap with O-ring for top compartment in HI 9033 housing

HI 740172 Protective cap for *Piccolo*

HI 740173 Protective cap for *Skincheck*™ electrode

HI 740178 Protective cap for RH probes

HI 740194 Protective cap for pH and ORP electrodes

HI 740208P Protective cap for FC 200 (2 pcs)

HI 740210P Protective cap for HI 98107 housing (2 pcs)

HI 740211 Transparent protective cap for D.O. probes

### **Cap, Protective, for Photometers**

Protects measurement cell of photometers from dust, liquid or any other damaging material, when not in use.

HI 740095 Protective cap for photometer measurement cell (2 pcs)



## Accessories



### Cap for Refillable Electrodes

Plastic cap for refilling hole of all pH and ORP electrodes using liquid electrolyte.  
HI 740093P Cap for refillable electrodes (20 pcs)

### Cap, Sintered, for RH Probes

Sintered cap for protecting the sensors of HANNA instruments® RH probes in harsh environments.  
HI 710011 Sintered cap

### Cap, Transportation, for D.O. Probes

Protects the membrane of HANNA instruments® D.O. probes during transportation, against accidental impacts.  
HI 740094P Transportation cap for D.O. probes (10 pcs)

### Carbon Dioxide, Test Kit Reagents

see Section V

### Carrying Case, Rugged

Rugged carrying case, 381x311x108 mm (15x12.2x4.2''), for photometers and turbidity meters with pre-cut sponge to fit instrument, reagents and apparatus.

HI 740218 Carrying case for portable photometers and turbidity meters

HI 740222 Carrying case for HI 832xx series

Rugged carrying case, 340x230x90 mm (13.4x9.1x3.5''), for general use.

HI 710031 General rugged carrying case

Rugged carrying case, 340x230x90 mm (13.4x9.1x3.5''), with thermoformed tray for HI 8014 and HI 8314 meters, complete with calibration screwdriver and pH 4 & pH 7 buffer solutions.

HI 721311 Kit for HI 8014 and HI 8314

Rugged carrying case, 340x230x90 mm (13.4x9.1x3.5''), with thermoformed tray for HI 931000 and all HANNA instruments® pH meters with HI 8424 casing, complete with pH 4 and pH 7 buffer solutions.

HI 721312 Kit for HI 991001 and HI 8424

Rugged carrying case, 340x230x90 mm (13.4x9.1x3.5''), with thermoformed tray for HI 8733 and HI 8633, complete with calibration screwdriver and HI 7030M calibration solution.

HI 721313 Kit for HI 8733 and HI 8633

Rugged carrying case, 340x230x90 mm (13.4x9.1x3.5''), with thermoformed tray for HI 8734, complete with calibration screwdriver and HI 7032M calibration solution.

HI 721314 Kit for HI 8734

Rugged carrying case, 340x230x90 mm (13.4x9.1x3.5''), with thermoformed tray for HI 8424 and HI 93510 casings.

HI 721316 Carrying case for portable meters

Rugged carrying case, 340x230x90 mm (13.4x9.1x3.5''), with thermoformed tray for HI 9025 casing and portable meters with built-in printer.

HI 721317 Carrying case for portable, waterproof and printing meters



## Accessories

Rugged carrying case, 340x230x90 mm (13.4x9.1x3.5"), for **HI 93703** turbidity meter, complete with calibration and cleaning solutions, tissues and cuvetts.  
**HI 731313** Maintenance kit for **HI 93703**

Rugged carrying case, 340x230x90 mm (13.4x9.1x3.5"), for **HI 93102** and **HI 93114**, complete with calibration and cleaning solutions, tissues and cuvetts.  
**HI 731327** Maintenance kit for **HI 93102** and **HI 93114**

Rugged carrying case with thermoformed tray for chemical test kits. The tray allows orderly placement of the reagents in each test kit.

**HI 740098** Carrying case for **HI 3813**

**HI 740099** Carrying case for **HI 3814**

**HI 740100** Carrying case for **HI 3816**

**HI 740101** Carrying case for **HI 3817**

**HI 740102** Carrying case for **HI 3819**

**HI 740103** Carrying case for **HI 3821**

**HI 740104** Carrying case for **HI 3823**

**HI 740105** Carrying case for **HI 3825**

**HI 740106** Carrying case for **HI 3827**

**HI 740107** Carrying case for **HI 3828**

**HI 740108** Carrying case for **HI 3837**

**HI 740115** Carrying case for **HI 3835**

**HI 740197** Carrying case for **HI 3810**

**HI 740201** Carrying case for **HI 3822**

In addition, several rugged carrying cases are available for other **HANNA** instruments' products, such as pocket meters.

**HI 740109** Carrying case for **Piccolo**

**HI 740110** Carrying case for **HI 98107** housing

**HI 740112** Carrying case for **Conmet**

### Carrying Case, Soft

Soft carrying case, 230x130x50 mm (9.1x5.1x2"), for **HI 8014** and portable photometer housing.

**HI 710001** Soft carrying case for **HI 8014** housing

Soft carrying case, 280x120x50 mm (11x4.7x2"), for **HI 9040** type thermistor thermometers.

**HI 710002** Soft carrying case for portable thermistor thermometers

Soft carrying case, 230x100x50 mm (9.1x3.9x2"), for **HI 9214** stick pH meter and any pocket meters, such as **Piccolo** and **Conmet**.

**HI 710003** Soft carrying case for **HI 9214**

Soft carrying case, 150x90x30 mm (5.9x3.5x1.2"), for **HI 93510** and **HI 95501** housings

**HI 710004** Soft carrying case for **HI 93510** housing

Soft carrying case for **HYGROCHECK®**.

**HI 740203** Soft carrying case for **HYGROCHECK®**

