

Accessories



- HI 740121 Right plastic cover for wall mounted instruments
- HI 740122P Plastic cover for terminals of **BL 7916** and **BL 7917** (2 pcs)
- HI 740124P Rear cover for **HI 140** series (2 pcs)
- HI 740125P Rear cover for **HI 981401** (2 pcs)
- HI 740126 Rear cover for **HI 8614** housing
- HI 740127 Transparent cover for built-in printer of portable meters
- HI 740129 Transparent front cover for **BL 7916** and **BL 7917**
- HI 740130 Transparent front cover for **HI 8710** housing

Cube, Color for CTK

Spare comparison cubes for colorimetric chemical test kits.

- HI 3824-99 Color cube for ammonia test kit
- HI 3830-99 Color cube for bromine test kit
- HI 3831-99 Color cube for chlorine test kit
- HI 3832-99 Color cube for iodine test kit
- HI 3833-99 Color cube for phosphate test kit
- HI 3834-99 Color cube for iron test kit
- HI 3874-99 Color cube for nitrate test kit

Cuvet

The glass cuvet is the most sensitive part for taking accurate measurements with photometers and turbidity meters. It is important that the cuvet is always clean and without scratches.

- HI 731321 Glass cuvet for **HI 937xx** portable meters & **HI 832xx** series (4 pcs)
- HI 731331 Glass cuvet for **HI 957xx** series (4 pcs)
- HI 731311 Glass vial with cap for **HI 93214** (5 pcs)

Cyanide, Colorimeter Reagents

see Section V

Cyanide, Test Kit Reagents

see Section V

Cyanuric Acid, Colorimeter Reagents

see Section V

Cyanuric Acid, Test Kit Reagents

see Section V

D

Demineralized Water

- Demi-02 Bottle to prepare 2 liters of demineralized water
- Demi-05 Bottle to prepare 5 liters of demineralized water
- Demi-10 Bottle to prepare 10 liters of demineralized water

Detergents, Colorimeter Reagents

see Section V

Detergents, Test Kit Reagents

see Section V

