

### E

#### Electrode Holder

Cantilever stand with a large steel base provides greater stability and flexibility, and holds up to 3 electrodes in any position.

- HI 76404 Electrode holder for pH 209 type bench meters
- HI 76405 Electrode holder with steel base
- HI 76405A Cantilever arm
- HI 76405C Steel base for HI 76405

Industrial electrode holder for immersion installation. Made of rugged PVC, this item can withstand most chemicals.

- HI 6050 Electrode holder for immersion installation, 0.6 m (23.8")
- HI 6051 Electrode holder for immersion installation, 1 m (3.3')
- HI 6052 Electrode holder for immersion installation, 1.6 m (5.25')

Industrial electrode holder for in-line installation. Two models are available, depending on the electrode thread. Made of fiber-glass reinforced polypropylene, these holders perfectly fit the relevant HANNA instruments® pH and ORP electrodes.

- HI 6054B Electrode holder for in-line installation, 3/4x16UNF thread
- HI 6054T Electrode holder for in-line installation, PG 13.5 thread

#### Electrode, ORP

see Section E

#### Electrode, pH

see Section E

#### Electrodes for Testers, Spare Parts

- HI 73106 Spare electrode for HI 98106 and HI 98107
- HI 73108 Spare electrode for HI 98108
- HI 73120 Spare electrode for HI 98120
- HI 73127 Spare electrode for HI 98121, HI 98127, HI 98128, HI 98129 and HI 98130
- HI 73201 Spare electrode for HI 98201
- HI 73202 Spare electrode for HI 98202 and HI 98203
- HI 1281 Spare electrode for Watercheck and Watercheck 1
- HI 1282 Spare electrode for Agritest and Agritest 1

#### Electrolyte Solution for D.O. Probes

see Solution, Electrolyte for D.O. Probes

#### Electrolyte Solution for pH and ORP Electrodes

see Solution, Electrolyte for pH and ORP Electrodes

### F

#### Filter for Chemical Test Kits

Replacement filter cartridge for HI 3838 formaldehyde test kit.

- HI 740135 Filter for HI 3838

