

## Accessories



### K

#### **Key for HI 141 Magnetic Start**

HI 740221 Key for HI 141 magnetic start

#### **Key for Testing Thermometer Calibration**

Allows a quick calibration check of the thermistor thermometers. Simply connect the key to the meter and verify that reading is within the key value  $\pm 0.4^{\circ}\text{C}$  (or  $\pm 0.8^{\circ}\text{F}$ ). If not, contact the nearest HANNA instruments® Customer Service Center for recalibration.

Test keys for thermometers using an HI 762 series probe.

HI 762-18C Test key at  $-18^{\circ}\text{C}$   
HI 762000C Test key at  $0^{\circ}\text{C}$   
HI 762070C Test key at  $70^{\circ}\text{C}$   
HI 762-004F Test key at  $-0.4^{\circ}\text{F}$   
HI 762032F Test key at  $32^{\circ}\text{F}$   
HI 762158F Test key at  $158^{\circ}\text{F}$

Test keys for thermometers using an HI 765 series probe.

HI 765-18C Test key at  $-18^{\circ}\text{C}$   
HI 765000C Test key at a  $0^{\circ}\text{C}$   
HI 765070C Test key at  $70^{\circ}\text{C}$   
HI 765-004F Test key at  $-0.4^{\circ}\text{F}$   
HI 765032F Test key at  $32^{\circ}\text{F}$   
HI 765158F Test key at  $158^{\circ}\text{F}$

### M

#### **Magnesium and Calcium, Test Kit Reagents**

see Section V

#### **Magnesium, Test Kit Reagents**

see Section V

#### **Manganese, Colorimeter Reagents**

see Section V

#### **Membrane for D.O. Probes**

This kit contains 5 ready-to-use, replacement Teflon® membranes.

HI 76407A/P Spare membrane for D.O. probes (5 pcs)

#### **Molybdenum, Colorimeter Reagents**

see Section V

#### **Motor for Magnetic Stirrers**

HI 740213 Motor for magnetic stirrers (3 pcs)

#### **Mounting Bracket**

see Bracket for Installation of Process Meters

