


## Accessories



### Power Cord

These cords power all **HANNA** instruments® meters requiring direct 115 Vac or 230 Vac supply.

- HI 740198 Power cord with DIN plug
- HI 740199 Power cord with US plug

### Power Plug

To be used with old wall mounted controllers and external devices (such as pumps or electrovalves) driven by the controllers, and with **BL 7916** & **BL 7917**.

- HI 767/P Input plug (5 pcs)
- HI 7671/P Output plug (5 pcs)

### Probe, Conductivity and TDS

#### Amperometric Probes

Spare graphite probes for EC and TDS testers.

- HI 73301 Spare probe for **HI 98301**
- HI 73302 Spare probe for **HI 98302**
- HI 73303 Spare probe for **HI 98303**
- HI 73304 Spare probe for **HI 98304**
- HI 73308 Spare probe for **HI 98308**
- HI 73309 Spare probe for **HI 98309**
- HI 73311 Spare probe for **HI 98311** and **HI 98312**

The following 2-pin probes have been designed for use with **Gro'chek** and **Pronto** series. They are amperometric probes, provided with built-in temperature sensor for automatic compensation, 2 m (6.6') cable and an Ultem® protective sleeve, which resists most aggressive chemicals. Models are available, for tank and flow-thru mounting applications.

Conductivity probes for **HI 983304** and **HI 983305**.

- HI 7631/2 EC probe with 1/2" thread
- HI 7631/3 EC probe with 3/4" thread

EC/TDS probes for **HI 983302**, **HI 983307** and **HI 983308**.

- HI 7630 EC/TDS probe
- HI 7632 EC/TDS probe
- HI 7632/2 EC/TDS probe with 1/2" thread
- HI 7632/3 EC/TDS probe with 3/4" thread

EC/TDS probes for **HI 983301**, **HI 983303** and **HI 983306**.

- HI 7634 EC/TDS probe
- HI 7634/2 EC/TDS probe with 1/2" thread
- HI 7634/3 EC/TDS probe with 3/4" thread

EC/TDS probe for **BL 983317**, **BL 983318** and **BL 983327** mini controllers.

- HI 7632-00 EC/TDS probe for high range mini controllers
- EC/TDS probe for **BL 983313**, **BL 983315**, **BL 983319**, **BL 983320**, **BL 983321**, **BL 983322**, **BL 983324** and **BL 983329** mini controllers.
- HI 7634-00 EC/TDS probe for low range mini controllers

EC/TDS probe for **HI 981409** and **HI 981410**.

- HI 1220 EC/TDS probe for fertigation stick meters

EC/TDS probe for **HI 9914**

- HI 7632D EC/TDS probe with DIN connector

