

S

Salinity Solution

see Sodium Chloride Solution

Salinity, Test Kit Reagents

see Section V

SAR, Test Kit Reagents

see Section V

Screwdriver for Calibration

Small calibration screwdriver, length 90 mm (3.5"), to perform manual calibration with trimmers.

HI 731326 Calibration screwdriver (20 pcs)

Silica, Colorimeter Reagents

see Section V

Silica, Test Kit Reagents

see Section V

Silver, Colorimeter Reagents

see Section V

Sleeve, Protective, for Conductivity Probes

These sleeves are needed to correctly perform EC and TDS measurements with HANNA instruments® probes. In addition, they protect sensors against accidental impact.

- HI 740179 Protective PVC sleeve for Conmet
- HI 740180 Protective PVC sleeve for HI 76302W
- HI 740181 Protective PVC sleeve for HI 1285
- HI 740182 Protective PVC sleeve for HI 76303
- HI 740189 Protective sleeve for HI 7685, HI 7686 and HI 7687

Sodium Chloride Solution

To be used for NaCl calibration of the relevant HANNA instruments® meters

- HI 7037L NaCl standard solution, 500 mL bottle
- HI 7037M NaCl standard solution, 230 mL bottle

Sodium Solution

see Section F

Software

Use the following Windows® compatible softwares to handle PC connection of the relevant HANNA instruments® meters.

- HI 92000 Windows® compatible software for logging instruments, photometers and turbidity meters
- HI 92140 Windows® compatible software for HI 140 and HI 142 series
- HI 92500 Windows® compatible software for microprocessor-based process controllers
- HI 800104 Windows® compatible software for HI 8000 controllers

