

Accessories



HI 981500 Windows® compatible software for pH Turtle (HI 9815)
HI 141000 Windows® compatible software for HI 141 series

Solution, Cal-Check™ for Photometers

see Section V

Solution for Cleaning Cuvets

For accurate measurements with photometers and turbidity meters, clean the glass cuvette before inserting it in the cell.

HI 93703-50 Cleaning solution for cuvettes, 230 mL bottle

Solution for Cleaning Electrodes

see Section F

Sulfate, Colorimeter Reagents

see Section V

Solution for EC and TDS Calibration

see Section F

Solution for Electrode Maintenance

see Section F

Solution for Electrode Storage

see Section F

Solution, Electrolyte for D.O. Probes

The regular maintenance of the D.O. probe is crucial to have accurate D.O. measurements.

HI 7041S Refilling electrolyte solution for D.O. probes (30 mL)

Solution, Electrolyte for pH and ORP Electrodes

see Section F

Solution for ORP

see Section F

Solution for pH Calibration

see Section F

Solution for Turbidity Calibration

For performing periodic calibration of turbidity meters, and ensure maximum accuracy. HANNA instruments® turbidity solutions do not contain unstable and carcinogenic formazine.

Solutions for HI 93102 and HI 93114

HI 93102-0 AMCO-AEPA turbidity solution at 0 NTU (30 mL)

HI 93102-20 AMCO-AEPA turbidity solution at 20 NTU (30 mL)

Solutions for HI 93703, HI 93703-11, LP 2000, LP 2000-11 and HI 93125

HI 93703-0 AMCO-AEPA turbidity solution at 0 FTU (30 mL)

HI 93703-10 AMCO-AEPA turbidity solution at 10 FTU (30 mL)

HI 93703-05 AMCO-EPA turbidity solution at 500 FTU (30 mL)

