

Spoon

Plastic spoon for precise dosing of powder chemical test kit reagents.

HI 740132P Plastic dosing spoon (10 pcs)

HI 740133 Plastic dosing spoon (6 pcs)

Sulfate, Colorimeter Reagents

see Section V

Sulfate, Test Kit Reagents

see Section V

Sulfite, Test Kit Reagents

see Section V

Syringe

Graduated syringe for accurately reagent dosing for CTK and photometers.

HI 740142P Graduated syringe (10 pcs)

HI 740143 Graduated syringe (6 pcs)

T

Temperature Probe

see Section P

Terminals, Extractable, for Process Instruments

These extractable modules have been designed for quick and easy wiring of HANNA instruments' process meters. They are available for 2, 3 and 4 inputs.

HI 740191P Extractable terminal, 2 inputs (10 pcs)

HI 740192P Extractable terminal, 3 inputs (10 pcs)

HI 740193P Extractable terminal, 4 inputs (10 pcs)

HI 740204P Large extractable terminal, 2 inputs (10 pcs)

HI 740205P Large extractable terminal, 3 inputs (10 pcs)

Test Tube

HI 740220 25 mL glass test tube (2 pcs)

Tissue for Wiping Cuvets

These non-abrasive tissues are ideal to wipe measurement cuvetts of photometers and turbidity meters, without damaging the cuvette surface.

HI 731318 Tissue for wiping cuvetts (4 pcs)

Titration System

Accessories for dosing reagents. The glass burette allows dosing of reagents with 0.1 mL resolution, and is supplied with a reagent glass container to perform up to 100 tests.

HI 740195 10 mL glass burette with 500 mL reagent container and pressure bulb

HI 740212 25 mL burette with 1 L reagent container and pressure bulb

HANNA instruments® provides a wide range of titration solutions.

HI 70455 0.01 N NaOH solution (1 L)

