

Reagents for Chemical Test Kits

Code	Description
Soil analysis, Ned (N), ascorbic acid (P), tetraphenylborate (K), pH indicators (pH)	
Range:	Nitrogen (N): trace, low, medium, high Phosphorus (P): trace, low, medium, high Potassium (K): trace, low, medium, high pH: 4 to 9 pH
Reagents	
HI 3895-010	Reagent kit for 10 tests (N, P, K, pH)
HI 3896-025	Reagent kit for 25 tests (N, P, K, pH)
Sulphate, barium chloride method	
Range:	20 to 100 mg/L
Reagents	
HI 38000-10	Reagent kit for 100 tests (SO ₄ ²⁻)
Sulphate LR/HR, barium chloride method	
Range:	100 to 1000 mg/L 1000 to 10000 mg/L
Reagents	
HI 38001-10	Reagent kit for 200 tests (SO ₄ ²⁻)
Sulphite, iodometric titration method	
Range:	0 to 20 mg/L 0 to 200 mg/L
Reagents	
HI 3822-100	Reagent kit for 100 tests (Na ₂ SO ₃)
Zinc, extended range, zincon method	
Range:	0.0 to 4.0 mg/L 0.0 to 20.0 mg/L
Reagents	
HI 38076-100	Reagent kit for 100 tests (Zn)
Zinc, zincon method	
Range:	0.0 to 3.0 mg/L
Reagents	
HI 3854-100	Reagent kit for 100 tests (Zn)

Accessories for Analysis

Accessories	
HI 3824-99	Ammonia color cube
HI 3830-99	Bromine color cube
HI 3831-99	Chlorine color cube
HI 3832-99	Iodine color cube
HI 3833-99	Phosphate color cube
HI 3834-99	Iron color cube
HI 731340	200 µL automatic pipet
HI 731341	1000 µL automatic pipet
HI 731342	2000 µL automatic pipet
HI 731350	Tip for 200 µL automatic pipet (25 pcs)
HI 731351	Tip for 1000 µL automatic pipet (25 pcs)
HI 731352	Tip for 2000 µL automatic pipet (25 pcs)
HI 740032P	Cap for 20 mL beaker (10 pcs)
HI 740034P	Cap for 50 mL beaker (10 pcs)
HI 740036P	Plastic Beaker, 50 mL (10 pcs)
HI 740037P	Plastic Beaker, 20 mL (10 pcs)
HI 740041	60 mL glass BOD bottle
HI 740048P	Plastic bottle, 120 mL (10 pcs)
HI 740049P	Cap for 120 mL plastic bottle (10 pcs)
HI 740132P	Dosing spoon (10 pcs)
HI 740142P	Graduated syringe (10 pcs)
HI 740144P	Tip for graduated syringe (10 pcs)
HI 740196P	Plastic vials (5 pcs)
HI 740220	25 mL glass test tube with cap, (2 pcs)
Demi-02	Bottle for 2 L of demineralized water
Demi-05	Bottle for 5 L of demineralized water
Demi-10	Bottle for 10 L of demineralized water

Demi-02, Demi-05 and Demi-10 Bottles for Demineralized Water

