

## Reagents for Turbidity Meters

Code	Description
<b>Bromine, method DPD</b>	
Range:	0.00 to 8.00 mg/L
Reagents	
HI 93716-01	Reagent kit for 100 tests (Br <sub>2</sub> )
HI 93716-03	Reagent kit for 300 tests (Br <sub>2</sub> )
<b>Chlorine, free, DPD method</b>	
Range:	0.00 to 2.50 mg/L
Reagents	
HI 93701-01	Powder reagent kit for 100 tests (free Cl <sub>2</sub> )
HI 93701-03	Powder reagent kit for 300 tests (free Cl <sub>2</sub> )
HI 93701-F	Liquid reagent kit for 300 tests (free Cl <sub>2</sub> )
<b>Chlorine, total, DPD method</b>	
Range:	0.00 to 3.50 mg/L
Reagents	
HI 93711-01	Powder reagent kit for 100 tests (total Cl <sub>2</sub> )
HI 93711-03	Powder reagent kit for 300 tests (total Cl <sub>2</sub> )
HI 93701-T	Liquid reagent kit for 300 tests (total Cl <sub>2</sub> )
HI 93711-D3	DPD3 reagent for 600 tests (total Cl <sub>2</sub> )
<b>Cyanuric Acid, turbidimetric method</b>	
Range:	0 to 80 mg/L
HI 93722-01	Reagent kit for 100 tests (cyanuric acid)
HI 93722-03	Reagent kit for 300 tests (cyanuric acid)
<b>Iodine, method DPD</b>	
Range:	0.0 to 12.5 mg/L
Reagents	
HI 93718-01	Reagent kit for 100 tests (I <sub>2</sub> )
HI 93718-03	Reagent kit for 300 tests (I <sub>2</sub> )
<b>Iron LR, TPTZ method</b>	
Range:	0.00 to 1.00 mg/L
Reagents	
HI 93746-01	Reagent kit for 50 tests (Fe LR)
HI 93746-03	Reagent kit for 150 tests (Fe LR)
<b>pH, phenol red method</b>	
Range:	5.9 to 8.5 pH
Reagents	
HI 93710-01	Reagent kit for 100 tests (pH)
HI 93710-03	Reagent kit for 300 tests (pH)
<b>Turbidity, ISO 7027 method</b>	
Range:	0.00 to 1000 FTU
Calibration Solutions	
HI 93703-0	AMCO-EPA solution @ 0 FTU (30 mL)
HI 93703-05	AMCO-EPA solution @ 500 FTU (30 mL)
HI 93703-10	AMCO-EPA solution @ 10 FTU (30 mL)
<b>Turbidity, USEPA 180.1 method</b>	
Range:	0.00 to 50.0 NTU
Calibration Solutions	
HI 93102-0	AMCO-EPA solution @ 0 NTU (30 mL)
HI 93102-20	SAMCO-EPA solution @ 20 NTU (30 mL)
HI 93124-03	calibration set @ 0.00, 2.50, 125 EBC (30 mL)
HI 93124-0	0.00 EBC solution (30 mL)
HI 93124-1	2.50 EBC solution (30 mL)
HI 93124-2	125 EBC solution (30 mL)

## Cal Check® Standard for HI 957xx Series

HI 95701-11	Cal Check™ Standard cuvetts for HI 95701
HI 95706-11	Cal Check™ Standard cuvetts for HI 95706
HI 95707-11	Cal Check™ Standard cuvetts for HI 95707
HI 95711-11	Cal Check™ Standard cuvetts for HI 95711 / HI 96711
HI 95715-11	Cal Check™ Standard cuvetts for HI 95715
HI 95717-11	Cal Check™ Standard cuvetts for HI 95717
HI 95721-11	Cal Check™ Standard cuvetts for HI 95721
HI 95724-11	Cal Check™ Standard cuvetts for HI 95724
HI 95729-11	Cal Check™ Standard cuvetts for HI 95729
HI 95734-11	Cal Check™ Standard cuvetts for HI 95734
HI 95747-11	Cal Check™ Standard cuvetts for HI 95747
HI 95761-11	Cal Check™ Standard cuvetts for HI 95761
HI 95762-11	Cal Check™ Standard cuvetts for HI 95762
HI 95769-11	Cal Check™ Standard cuvetts for HI 95769
HI 95771-11	Cal Check™ Standard cuvetts for HI 95771
HI 93701-00	Free Chlorine Standard 0 mg/L (4 pcs)
HI 93701-10	Free Chlorine Standard 1 mg/L (4 pcs)
HI 93701-11	Free Chlorine Standard 0, 0.25, 0.5 and 1 mg/L (1 pc each)
HI 93701-25	Free Chlorine Standard 0.25 mg/L (4 pcs)
HI 93701-50	Free Chlorine Standard 0.5 mg/L (4 pcs)



**Cal Check™ Standard Solutions:**  
**Calibrate and Check HI 957xx series!**